

**510(k) Summary**  
**Terrats Medical SL**  
**DESS® Dental Smart Solutions**  
**K212577**

**December 10, 2021**

**ADMINISTRATIVE INFORMATION**

Manufacturer Name	Terrats Medical SL Carrer Mogoda, 75-99 08210 Barberà del Vallès Barcelona, Spain
Telephone	+34 935 646 006
Official Contact	Roger Terrats, CEO
Representative/Consultant	Melissa Burbage, BS Floyd G. Larson, MS, MBA Kevin A. Thomas, PhD PaxMed International, LLC 12264 El Camino Real, Suite 400 San Diego, CA 92130 Telephone +1 858-792-1235 Fax +1 858-792-1236 Email mburbage@paxmed.com flarson@paxmed.com kthomas@paxmed.com

**DEVICE NAME AND CLASSIFICATION**

Trade/Proprietary Name	DESS Dental Smart Solutions
Common Name	Dental implant abutment
Regulation Number	21 CFR 872.3630
Regulation Name	Endosseous dental implant abutment
Regulatory Class	Class II
Product Code	NHA
Classification Panel	Dental Products Panel
Reviewing Division	DHT1B: Division of Dental Devices

**PREDICATE DEVICE INFORMATION**

Primary Predicate Device  
K191986, DESS Dental Smart Solutions, Terrats Medical SL

Additional Predicate Devices  
K170588, DESS Dental Smart Solutions, Terrats Medical SL  
K173908, DESS Dental Smart Solutions, Terrats Medical SL

Reference Devices  
K193425, Pre-Milled Blank, ARUM Dentistry Co. Ltd  
K183518, Preat Abutments, Preat Corporation

## INDICATIONS FOR USE STATEMENT

DESS Dental Smart Solutions abutments are intended to be used in conjunction with endosseous dental implants in the maxillary or mandibular arch to provide support for prosthetic restorations.

All digitally designed custom abutments for use with Pre-milled Blanks are to be sent to a Terrats Medical validated milling center for manufacture.

## Compatible Implant Systems

Compatible Implant System	Implant Body Diameter, mm	Implant Platform
NobelActive®, NobelParallel Conical	3.5	NP
	4.3, 5.0	RP
	5.5	WP
Straumann® Bone Level	3.3	NC
	4.1/4.8	RC
Zimmer Screw-Vent®/ Tapered Screw-Vent®	3.7, 4.1	3.5
	4.7	4.5
	6.0	5.7

## SUBJECT DEVICE DESCRIPTION

The purpose of this submission is to implement a labeling change for certain Pre-milled Blanks of the DESS Dental Smart Solutions abutment system, cleared under K170588, K173908 and K191986, to permit the manufacture of custom abutments with angulation up to 30° and to add three (3) Pre-milled Blanks that also will include such labeling. Design parameters cleared in K170588, K173908 and K191986 were limited to straight abutments only. This change in labeling is for Pre-milled Blanks compatible with three (3) systems: NobelActive/NobelParallel Conical Connection, Straumann Bone Level, and Zimmer Screw-Vent/Tapered Screw-Vent implants. Note that, because NobelActive and Nobel Parallel Conical Connection share the same implant/abutment interface, they are considered one system for purposes of this submission, as they were in K170588, K173908 and K191986. No new compatibilities are added in this submission.

Pre-milled Blanks are designed for fabrication of a custom abutment by a CAD/CAM process. All patient-specific custom abutment fabrication is by prescription on the order of the clinician. Nobel-compatible Pre-milled Blanks are made of titanium alloy (Ti-6Al-4V) conforming to ASTM F136 *Standard Specification for Wrought Titanium-6Aluminum-4Vanadium ELI (Extra Low Interstitial) Alloy for Surgical Implant Applications (UNS R56401)* or from cobalt chromium alloy (Co-Cr-Mo) conforming to ASTM F1537 *Standard Specification for Wrought Cobalt-28Chromium-6Molybdenum Alloys for Surgical Implants (UNS R31537, UNS R31538, and UNS R31539)*. Straumann-compatible and Zimmer-compatible Pre-milled Blanks are made of titanium alloy (Ti-6Al-4V). They are available in engaging designs and are compatible with the implant systems shown the table below *Compatible Implant Systems and Platforms*. Except for the maximum angulation of the final abutment and the addition of three (3) Pre-milled Blanks of 14 mm diameter to the previously cleared Pre-milled Blanks of 10 mm and 14 mm diameter, subject device Pre-milled Blanks are identical to Pre-milled Blanks cleared (with slight variations on the name) in K170588, K173908 and K191986. Compatibility with the implant platform was demonstrated in K170588 and K191986.

The design parameters for the CAD/CAM fabrication of custom abutments from subject device Pre-milled Blanks are:

- Minimum wall thickness – 0.45 mm
- Minimum post height – 4.0 mm

Maximum gingival height – 6.0 mm  
 Minimum gingival height – 0.3 mm  
 Maximum angulation of the final abutment – 30°

Manufacture of CAD/CAM custom abutments from Pre-milled Blanks is to be performed at a Terrats Medical validated milling center.

#### Compatible Implant Systems and Platforms

Compatible Implant System	DESS Abutment System	Implant Body Diameter, mm	Implant Platform	Connection
NobelActive®, NobelParallel Conical	Active Hex	3.5	NP	Internal
		4.3, 5.0	RP	
		5.5	WP	
Straumann® Bone Level	Conical BL	3.3	NC	Internal
		4.1/4.8	RC	
Zimmer Screw-Vent®/ Tapered Screw-Vent®	Internal Hex USA	3.7, 4.1	3.5	Internal
		4.7	4.5	
		6.0	5.7	

#### PERFORMANCE DATA

Non-clinical data submitted to demonstrate substantial equivalence included: static and dynamic testing according to ISO 14801. Information was leveraged from prior clearance (K170588 and K173908) to demonstrate substantial equivalence with regards to sterilization and biocompatibility. No clinical data were included in this submission.

#### EQUIVALENCE TO MARKETED DEVICES

The subject device is substantially equivalent in indications and design principles to the primary predicate device and the additional predicate devices listed above. Provided at the end of this summary are tables comparing the Indications for Use Statements and the technological characteristics of the subject device, the primary predicate device, the additional predicate devices, and the reference devices.

Subject device abutments are substantially equivalent in intended use to the primary predicate device cleared in K191986, the additional predicate devices K170588, K173908 and the reference device K193425. All are intended for use with endosseous dental implants to provide functional and esthetic rehabilitation of the edentulous maxilla and mandible. The Indications for Use Statement (IFUS) for the subject device is substantially equivalent to that of the primary predicate K191986 and additional predicate devices K170588 and K173908.

All subject device abutments are identical in design, materials, and technological characteristics to those of the primary predicate K191986 and additional predicate devices K170588 and K173908. The only change to the subject device is a labeling change to permit custom abutments to be manufactured from the Pre-milled Blanks with an angulation up to 30°, which is substantially equivalent to that of the reference device K193425, and the addition of three new Pre-milled Blanks of 14 mm diameter to the Pre-milled Blanks of 10 mm and 14 mm diameter previously cleared for the predicate devices. No new compatibilities are included in this submission.

Substantial equivalence of the subject device with design parameters that permit angulation of the final abutment is supported by dynamic testing according to ISO 14801 *Dentistry – Implants – Dynamic fatigue test for endosseous dental implants* and by the reference device K183518.

The subject device is to be sterilized by the end-user, the same as primary predicate device K191986 and additional predicate devices K170588 and K173908.

**CONCLUSION**

The subject device, the primary predicate device, additional predicate devices, and the reference devices have the same intended use, have similar technological characteristics, and are made of the same materials. The subject device, the primary predicate, and the additional predicate devices encompass the same range of physical dimensions, are packaged in similar materials, and are to be sterilized using similar methods. The data included in this submission demonstrate substantial equivalence to the predicate devices listed above.

**Table of Substantial Equivalence – Indications for Use Statement**

<p align="center"><b>Subject Device</b></p> <p align="center"><b>DESS Dental Smart Solutions</b></p> <p align="center"><b>K212577</b></p> <p align="center"><b>Terrats Medical SL</b></p>	<p>DESS Dental Smart Solutions abutments are intended to be used in conjunction with endosseous dental implants in the maxillary or mandibular arch to provide support for prosthetic restorations.</p> <p>All digitally designed custom abutments for use with Pre-milled Blanks are to be sent to a Terrats Medical validated milling center for manufacture.</p> <p align="center">Compatible Implant Systems</p> <table border="1"> <thead> <tr> <th>Compatible Implant System</th> <th>Implant Body Diameter, mm</th> <th>Implant Platform</th> </tr> </thead> <tbody> <tr> <td rowspan="3">NobelActive®, NobelParallel Conical</td> <td>3.5</td> <td>NP</td> </tr> <tr> <td>4.3, 5.0</td> <td>RP</td> </tr> <tr> <td>5.5</td> <td>WP</td> </tr> <tr> <td rowspan="2">Straumann® Bone Level</td> <td>3.3</td> <td>NC</td> </tr> <tr> <td>4.1/4.8</td> <td>RC</td> </tr> <tr> <td rowspan="3">Zimmer Screw Vent®/ Tapered Screw-Vent®</td> <td>3.7, 4.1</td> <td>3.5</td> </tr> <tr> <td>4.7</td> <td>4.5</td> </tr> <tr> <td>6.0</td> <td>5.7</td> </tr> </tbody> </table>	Compatible Implant System	Implant Body Diameter, mm	Implant Platform	NobelActive®, NobelParallel Conical	3.5	NP	4.3, 5.0	RP	5.5	WP	Straumann® Bone Level	3.3	NC	4.1/4.8	RC	Zimmer Screw Vent®/ Tapered Screw-Vent®	3.7, 4.1	3.5	4.7	4.5	6.0	5.7																												
Compatible Implant System	Implant Body Diameter, mm	Implant Platform																																																	
NobelActive®, NobelParallel Conical	3.5	NP																																																	
	4.3, 5.0	RP																																																	
	5.5	WP																																																	
Straumann® Bone Level	3.3	NC																																																	
	4.1/4.8	RC																																																	
Zimmer Screw Vent®/ Tapered Screw-Vent®	3.7, 4.1	3.5																																																	
	4.7	4.5																																																	
	6.0	5.7																																																	
<p align="center"><b>Primary Predicate Device</b></p> <p align="center"><b>K191986</b></p> <p align="center"><b>DESS Dental Smart Solutions</b></p> <p align="center"><b>Terrats Medical SL</b></p>	<p>DESS Dental Smart Solutions abutments are intended to be used in conjunction with endosseous dental implants in the maxillary or mandibular arch to provide support for prosthetic restorations.</p> <p>All digitally designed custom abutments for use with Ti Base abutments or Pre-milled (Blank) abutments are to be sent to a Terrats Medical validated milling center for manufacture.</p> <p align="center">Compatible Implant Systems</p> <table border="1"> <thead> <tr> <th>Compatible Implant System</th> <th>Implant Body Diameter, mm</th> <th>Implant Platform</th> </tr> </thead> <tbody> <tr> <td>Ankylos C/X</td> <td>3.5, 4.5, 5.5</td> <td>2.52 mm</td> </tr> <tr> <td rowspan="3">Astra Tech EV</td> <td>3.6</td> <td>2.9 mm</td> </tr> <tr> <td>4.2</td> <td>3.5 mm</td> </tr> <tr> <td>4.8</td> <td>4.1 mm</td> </tr> <tr> <td rowspan="3">Astra Tech OsseoSpeed™</td> <td>3.0</td> <td>3.0 mm</td> </tr> <tr> <td>3.5/4.0</td> <td>3.5/4.0 mm</td> </tr> <tr> <td>4.5/5.0</td> <td>4.5/5.0 mm</td> </tr> <tr> <td rowspan="3">Biomet 3i Certain®</td> <td>3.25</td> <td>3.45 mm</td> </tr> <tr> <td>4.0</td> <td>4.1 mm</td> </tr> <tr> <td>5.0</td> <td>5.0 mm</td> </tr> <tr> <td rowspan="3">Biomet 3i OSSEOTITE®</td> <td>3.25</td> <td>3.4 mm</td> </tr> <tr> <td>3.75, 4.0</td> <td>4.1 mm</td> </tr> <tr> <td>5.0</td> <td>5.0 mm</td> </tr> <tr> <td rowspan="4">Camlog</td> <td>3.3</td> <td>3.3 mm</td> </tr> <tr> <td>3.8</td> <td>3.8 mm</td> </tr> <tr> <td>4.3</td> <td>4.3 mm</td> </tr> <tr> <td>5.0</td> <td>5.0 mm</td> </tr> <tr> <td rowspan="3">FRIADENT XiVE®</td> <td>3.4</td> <td>3.4 mm</td> </tr> <tr> <td>3.8</td> <td>3.8 mm</td> </tr> <tr> <td>4.5</td> <td>4.5 mm</td> </tr> </tbody> </table>	Compatible Implant System	Implant Body Diameter, mm	Implant Platform	Ankylos C/X	3.5, 4.5, 5.5	2.52 mm	Astra Tech EV	3.6	2.9 mm	4.2	3.5 mm	4.8	4.1 mm	Astra Tech OsseoSpeed™	3.0	3.0 mm	3.5/4.0	3.5/4.0 mm	4.5/5.0	4.5/5.0 mm	Biomet 3i Certain®	3.25	3.45 mm	4.0	4.1 mm	5.0	5.0 mm	Biomet 3i OSSEOTITE®	3.25	3.4 mm	3.75, 4.0	4.1 mm	5.0	5.0 mm	Camlog	3.3	3.3 mm	3.8	3.8 mm	4.3	4.3 mm	5.0	5.0 mm	FRIADENT XiVE®	3.4	3.4 mm	3.8	3.8 mm	4.5	4.5 mm
Compatible Implant System	Implant Body Diameter, mm	Implant Platform																																																	
Ankylos C/X	3.5, 4.5, 5.5	2.52 mm																																																	
Astra Tech EV	3.6	2.9 mm																																																	
	4.2	3.5 mm																																																	
	4.8	4.1 mm																																																	
Astra Tech OsseoSpeed™	3.0	3.0 mm																																																	
	3.5/4.0	3.5/4.0 mm																																																	
	4.5/5.0	4.5/5.0 mm																																																	
Biomet 3i Certain®	3.25	3.45 mm																																																	
	4.0	4.1 mm																																																	
	5.0	5.0 mm																																																	
Biomet 3i OSSEOTITE®	3.25	3.4 mm																																																	
	3.75, 4.0	4.1 mm																																																	
	5.0	5.0 mm																																																	
Camlog	3.3	3.3 mm																																																	
	3.8	3.8 mm																																																	
	4.3	4.3 mm																																																	
	5.0	5.0 mm																																																	
FRIADENT XiVE®	3.4	3.4 mm																																																	
	3.8	3.8 mm																																																	
	4.5	4.5 mm																																																	

			5.5	5.5 mm
	MegaGen AnyRidge		3.5, 4.0, 4.5, 5.0, 5.5	3.5 mm
	NobelActive®, NobelReplace/NobelParallel Conical		3.0	3.0 (3.0 mm)
			3.5	NP (3.5 mm)
			4.3, 5.0	RP (3.9 mm)
			5.5	WP (5.1 mm)
	NobelReplace® Trilobe		3.5	NP (3.5 mm)
			4.3	RP (4.3 mm)
			5.0	WP (5.0 mm)
			6.0	6.0 (6.0 mm)
	Nobel Brånemark System®		3.3	NP (3.5 mm)
			3.75, 4.0	RP (4.1 mm)
			5.0	WP (5.1 mm)
	Osstem TS		3.5	Mini (2.8 mm)
			4.0, 4.5, 5.0, 6.0, 7.0	Regular (3.35 mm)
	Straumann® Bone Level		3.3	NC (3.3 mm)
			4.1/4.8	RC (4.1/4.8 mm)
	Straumann® Tissue Level		3.3	NNC (3.5 mm)
			3.3, 4.1, 4.8	RN (4.8 mm)
			4.8	WN (6.5 mm)
	Zimmer Screw Vent®/ Tapered Screw-Vent®		3.3, 3.7, 4.1	3.5 mm
			4.7	4.5 mm
			6.0	5.7 mm

<p><b>Additional Predicate Device</b></p> <p><b>K170588</b></p> <p><b>DESS Dental Smart Solutions</b> <b>Terrats Medical SL</b></p>	<p>DESS Dental Smart Solutions abutments are intended to be used in conjunction with endosseous dental implants in the maxillary or mandibular arch to provide support for prosthetic restorations.</p> <p>All digitally designed custom abutments for use with Ti Base or Pre-milled (Blank) abutments are to be sent to a Terrats Medical validated milling center for manufacture.</p> <p>Compatible Implant Systems</p>		
	Implant System Compatibility	Implant Diameter (mm)	Platform Diameter (mm)
	3i Certain®	3.25, 4.0, 5.0	3.4, 4.1, 5.0
	3i OSSEOTITE®	3.25, 3.75, 4.0, 5.0	3.4, 4.1, 5.0
	OsseoSpeed™	3.5, 4.0, 5.0	3.5/4.0, 4.5/5.0
	FRIADENT XiVE	3.4, 3.8, 4.5	3.4, 3.8, 4.5
	NobelActive®	3.5, 4.3, 5.0	NP, RP
	NobelReplace Conical	3.5, 4.3, 5.0	NP, RP
	Nobel Replace Trilobe	3.5, 4.3, 5.0	NP, RP, WP
	Brånemark	3.5, 3.75/4.0, 5.0	NP, RP, WP
	Straumann® Bone Level	3.3, 4.1, 4.8	NC, RC
	Straumann® Tissue Level	3.3, 4.1, 4.8	RN, WN
	Tapered Screw-Vent®	3.7, 4.1, 4.7, 6.0	3.5, 4.5, 5.7

<p style="text-align: center;"><b>Additional Predicate Device</b></p> <p style="text-align: center;"><b>K173908</b></p> <p style="text-align: center;"><b>DESS Dental Smart Solutions</b> <b>Terrats Medical SL</b></p>	<p>DESS Dental Smart Solutions abutments are intended to be used in conjunction with endosseous dental implants in the maxillary or mandibular arch to provide support for prosthetic restorations. All digitally designed custom abutments for use with Aurum™ Abutment or Pre-milled Blank are to be sent to a Terrats Medical validated milling center for manufacture.</p> <p>Compatible Implant Systems</p> <table border="1" data-bbox="772 344 1806 768"> <thead> <tr> <th>Implant System Compatibility</th> <th>Implant Body</th> <th>Implant Platform</th> </tr> </thead> <tbody> <tr> <td>3i Certain®</td> <td>3.25, 4.0, 5.0</td> <td>3.4, 4.1, 5.0</td> </tr> <tr> <td>3i OSSEOTITE®</td> <td>3.25, 3.75, 4.0, 5.0</td> <td>3.4, 4.1, 5.0</td> </tr> <tr> <td>OsseoSpeed™</td> <td>3.5, 4.0, 5.0</td> <td>3.5/4.0, 4.5/5.0</td> </tr> <tr> <td>FRIADENT XiVE</td> <td>3.4, 3.8, 4.5</td> <td>3.4, 3.8, 4.5</td> </tr> <tr> <td>NobelActive®</td> <td>3.5, 4.3, 5.0</td> <td>NP, RP</td> </tr> <tr> <td>NobelReplace® Conical</td> <td>3.5, 4.3, 5.0</td> <td>NP, RP</td> </tr> <tr> <td>NobelReplace® Trilobe</td> <td>3.5, 4.3, 5.0</td> <td>NP, RP, WP</td> </tr> <tr> <td>Brånemark</td> <td>3.5, 3.75/4.0, 5.0</td> <td>NP, RP, WP</td> </tr> <tr> <td>Straumann® Bone Level</td> <td>3.3, 4.1, 4.8</td> <td>NC, RC</td> </tr> <tr> <td>Straumann® Tissue Level</td> <td>3.3, 4.1, 4.8</td> <td>RP, WP</td> </tr> <tr> <td>Tapered Screw-Vent®</td> <td>3.7, 4.1, 4.7, 6.0</td> <td>3.5, 4.5, 5.7</td> </tr> </tbody> </table>	Implant System Compatibility	Implant Body	Implant Platform	3i Certain®	3.25, 4.0, 5.0	3.4, 4.1, 5.0	3i OSSEOTITE®	3.25, 3.75, 4.0, 5.0	3.4, 4.1, 5.0	OsseoSpeed™	3.5, 4.0, 5.0	3.5/4.0, 4.5/5.0	FRIADENT XiVE	3.4, 3.8, 4.5	3.4, 3.8, 4.5	NobelActive®	3.5, 4.3, 5.0	NP, RP	NobelReplace® Conical	3.5, 4.3, 5.0	NP, RP	NobelReplace® Trilobe	3.5, 4.3, 5.0	NP, RP, WP	Brånemark	3.5, 3.75/4.0, 5.0	NP, RP, WP	Straumann® Bone Level	3.3, 4.1, 4.8	NC, RC	Straumann® Tissue Level	3.3, 4.1, 4.8	RP, WP	Tapered Screw-Vent®	3.7, 4.1, 4.7, 6.0	3.5, 4.5, 5.7
Implant System Compatibility	Implant Body	Implant Platform																																			
3i Certain®	3.25, 4.0, 5.0	3.4, 4.1, 5.0																																			
3i OSSEOTITE®	3.25, 3.75, 4.0, 5.0	3.4, 4.1, 5.0																																			
OsseoSpeed™	3.5, 4.0, 5.0	3.5/4.0, 4.5/5.0																																			
FRIADENT XiVE	3.4, 3.8, 4.5	3.4, 3.8, 4.5																																			
NobelActive®	3.5, 4.3, 5.0	NP, RP																																			
NobelReplace® Conical	3.5, 4.3, 5.0	NP, RP																																			
NobelReplace® Trilobe	3.5, 4.3, 5.0	NP, RP, WP																																			
Brånemark	3.5, 3.75/4.0, 5.0	NP, RP, WP																																			
Straumann® Bone Level	3.3, 4.1, 4.8	NC, RC																																			
Straumann® Tissue Level	3.3, 4.1, 4.8	RP, WP																																			
Tapered Screw-Vent®	3.7, 4.1, 4.7, 6.0	3.5, 4.5, 5.7																																			
<p style="text-align: center;"><b>Reference Device</b></p> <p style="text-align: center;"><b>K193425</b></p> <p style="text-align: center;"><b>Pre-Milled Blank</b> <b>ARUM Dentistry Co. Ltd</b></p>	<p>ARUM Dentistry’s Pre-Milled Blank abutments are intended for attachment to dental implants in order to provide support for customized prosthetic restorations. Pre-Milled Blank abutments are indicated for screw-retained single restorations or cement-retained single or multi-unit restorations. The customized Pre-Milled Blank abutment will be attached to a dental implant using the included ARUM Dentistry prosthetic screw.</p> <p>ARUM Dentistry’s Pre-Milled Blanks are compatible with the implant systems listed in the Compatibility Table:</p> <p>Compatibility Table</p> <table border="1" data-bbox="697 1015 1879 1214"> <thead> <tr> <th colspan="2">ARUM Pre-Milled Blank</th> <th rowspan="2">Implant Platform compatibility</th> <th rowspan="2">Restorative Platform Diameter (mm)</th> <th rowspan="2">Implant Body diameter (mm)</th> <th rowspan="2">Abutment Screw</th> </tr> <tr> <th>10 mm</th> <th>14 mm</th> </tr> </thead> <tbody> <tr> <td>CIHE037</td> <td>CIHE038</td> <td>NobelActive NP</td> <td>3.5</td> <td>3.5</td> <td rowspan="2">CSTO001</td> </tr> <tr> <td>CIHE039</td> <td>CIHE040</td> <td>NobelActive NP</td> <td>3.9</td> <td>4.3/5.0</td> </tr> <tr> <td>CIHE135</td> <td>CIHE136</td> <td>NobelActive WP</td> <td>5.1</td> <td>5.5</td> <td>CSTO002</td> </tr> </tbody> </table> <p>All digitally-designed Pre-Milled Blank abutments are intended to be sent to an ARUM Dentistry-validated milling center for manufacture.</p>	ARUM Pre-Milled Blank		Implant Platform compatibility	Restorative Platform Diameter (mm)	Implant Body diameter (mm)	Abutment Screw	10 mm	14 mm	CIHE037	CIHE038	NobelActive NP	3.5	3.5	CSTO001	CIHE039	CIHE040	NobelActive NP	3.9	4.3/5.0	CIHE135	CIHE136	NobelActive WP	5.1	5.5	CSTO002											
ARUM Pre-Milled Blank		Implant Platform compatibility	Restorative Platform Diameter (mm)					Implant Body diameter (mm)	Abutment Screw																												
10 mm	14 mm																																				
CIHE037	CIHE038	NobelActive NP	3.5	3.5	CSTO001																																
CIHE039	CIHE040	NobelActive NP	3.9	4.3/5.0																																	
CIHE135	CIHE136	NobelActive WP	5.1	5.5	CSTO002																																

**Reference Device**  
**K183518**  
**Preat Abutments**  
**Preat Corporation**

Preat Abutments are intended to be used in conjunction with endosseous dental implants in the maxillary or mandibular arch to provide support for single-unit or multi-unit prosthetic restorations. The Titanium Base abutments consists of two major parts. Specifically, the titanium base and mesostructured components make up a two-piece abutment.  
All digitally designed custom abutments, superstructures, and/or hybrid crowns for use with Titanium Base or Titanium Blank are to be sent to a Preat validated milling center for manufacture.  
Compatible Implant Systems

Compatible Implant System	Implant Body Diameter (mm)	Implant Platform Diameter (mm)
3i OSSEOTITE® Certain®	3.25	3.4
	4.0	4.1
	5.0	5.0
	6.0	6.0
Astra Tech OsseoSpeed™	3.0	3.0
	3.5, 4.0	3.5/4.0
	4.5, 5.0	4.5/5.0
BioHorizons Tapered Internal	3.0	3.0
	3.5	3.5
	4.0	4.5
HIOSEN ET III	3.5	Mini
	4.0, 4.5, 5.0, 6.0, 7.0	Regular
Implant Direct Legacy	3.2	3.0
	3.7, 4.2	3.5
	4.7, 5.2	4.5
	5.7, 7.0	5.7
MegaGen AnyRidge	3.5, 4.0, 4.5, 5.0, 5.5	3.5
Neoss	3.5, 4.0, 4.5, 5.0, 5.5	4.1
NobelActive®	3.5	NP
	4.3, 5.0	RP
Nobel Replace™	3.5	NP
	4.0, 4.3, 5.0	RP
	5.0	WP
	6.0	6.0
Straumann® Bone Level	3.3	NC
	4.1, 4.8	RC
Straumann® Tissue Level	3.3, 4.1, 4.8	RN
	4.8, 6.5	WN
Zimmer Screw-Vent®/Tapered Screw-Vent®	3.3, 3.7, 4.1	3.5
	4.7	4.5
	6.0	5.7



**Table of Substantial Equivalence – Technological Characteristics**

Comparison	Subject Device	Primary Predicate Device	Additional Predicate Devices		Reference Devices	
	K212577 DESS Dental Smart Solutions Terrats Medical SL	K191986 DESS Dental Smart Solutions Terrats Medical SL	K170588 DESS Dental Smart Solutions Terrats Medical SL	K173908 DESS Dental Smart Solutions Terrats Medical SL	K193425 Pre-Milled Blank ARUM Dentistry Co. Ltd	K183518 Preat Abutments Preat Corporation
<b>Reason for Predicate</b>	n/a	Design	Design	Design	30° angle	Performance Testing
<b>Design</b>						
Abutment Designs	CAD/CAM Blanks	Healing, Temporary, Straight, Multi-unit, Locator-type, CAD/CAM Bases, CAD/CAM Blanks	Healing, Temporary, Straight, Multi-unit, Locator-type, CAD/CAM Bases, CAD/CAM Blanks	CAD/CAM Bases, CAD/CAM Blanks	CAD/CAM Blanks	Healing, Temporary, Straight, Multi-unit, CAD/CAM Bases, CAD/CAM Blanks
Prosthesis Attachment	Cement-retained Screw-retained	Cement-retained Screw-retained	Cement-retained Screw-retained	Cement-retained Screw-retained	Cement-retained	Cement-retained Screw-retained
Restoration	Single-unit, Multi-unit	Single-unit, Multi-unit	Single-unit, Multi-unit	Single-unit, Multi-unit	Single-unit, Multi-unit	Single-unit, Multi-unit
Abutment/Implant Platform Diameter, mm	2.52 – 6.0	2.52 – 6.0	3.0 – 6.0	3.6 - 5.0	3.5 – 5.5	3.0 – 6.5
Prosthetic Platform Diameter, mm	4.5 – 6.5	4.5-6.5	4.5	4.0 – 6.5	3.5 – 5.1	3.0 – 6.5
Blank Abutment Angulation	0° to 30°	0°	0°	0°	0° to 30°	0° to 30°
Abutment/ Implant Interface	Internal	Internal	Internal	Internal	Internal	Internal
<b>Material</b>	Ti-6Al-4V ELI Co-Cr-Mo Alloy	Ti-6Al-4V ELI Co-Cr-Mo Alloy	Ti-6Al-4V ELI	Ti-6Al-4V ELI Co-Cr-Mo Alloy	Ti-6Al-4V ELI	Ti-6Al-4V alloy